


SARS-CoV-2 Antibody Test (Lateral flow method)

- As an auxiliary and supplementary method for the further diagnosis to the patients with negative result in PCR.
- Easy to use, instant result in 15 minutes



According to the Novel Coronavirus Pneumonia Diagnosis and Treatment Plan (7th Edition) published by general office of national health committee, serology testing can be used as a diagnostic method of COVID-19 infection.



Antibodies can be detected as early as one week after infection, with slightly difference between individuals.

1. Negative result indicates no infection or infection at a very early stage (within first few days).
2. Positive result indicates an infection with coronavirus.

Order information

Product name	Packing size	Sample type	Storage condition	Shelf life
SARS-CoV-2Antibody Test (Lateral flow method)	20T	Whole blood, serum and plasma	2°C~30°C	6 months

