

Product Information	Product Description	Product Specifications				
Product Part Number	Enhanced Description	Label Size (WxH)	Colour	Core Size	Quantity/ unit	Cena bez DPH
LWCS-501-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	2" x 1" / 50.8mm x 25.4mm	White	1 inch	1000 labels/roll	3 012 Kč
LWCS-502-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	2" x 0.25" / 50.8mm x 6.4mm	White	1 inch	1000 labels/roll	2 496 Kč
LWCS-503-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	1.5" x 0.75" / 38.1mm x 19.1mm	White	1 inch	750 labels/roll	1 727 Kč
LWCS-504-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	1" x 1" / 25.4mm x 25.4mm	White	1 inch	750 labels/roll	1 727 Kč
LWCS-505-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	1" x 0.5" / 25.4mm x 12.7mm	White	1 inch	750 labels/roll	1 727 Kč
LWCS-506-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	0.94" x 0.5" + 0.375" / 23.9mm x 12.7mm + 9.5mm	White	1 inch	750 labels/roll	2 208 Kč
LWCS-507-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	1.125" x 0.625" + 0.437" / 28.6mm x 15.9mm + 11.1mm	White	1 inch	750 labels/roll	2 208 Kč
LWCS-510-WH	Our CryoSTUCK® labels for frozen surfaces are now available for printing with DYMO LabelWriter 5-Series printers. These durable, solvent-resistant cryogenic labels are designed for direct application to frozen surfaces as low as (-80°C/-112°F), including on vials and tubes stored in ultra-low temperature freezers. They offer a reliable way to label frozen containers without needing to thaw samples first, ensuring the integrity of valuable temperature-sensitive specimens.	0.5" / 12.7 mm (circle)	White	1 inch	2000 labels/roll	1 849 Kč

DYMO-LW550-EMEA	<p>The DYMO LabelWriter™ 550 is a direct thermal printer that does not require costly inks or toner to print high-quality text and barcodes, while its compact design allows it to fit seamlessly into any workspace. Its patented automatic label recognition technology will automatically identify label size and type, allowing it to quickly generate the template once you load the labels, so you never have to waste time selecting the right template. In addition, the automatic label counter indicates how many labels remain—so you'll never run out unexpectedly. Free DYMO Connect software is included with every printer, allowing you to create 1D and 2D barcodes along with crisp, clear text and graphics. The 550 is up to 20% faster than its predecessor, the LabelWriter™ 450, and is fully compatible with our LWCS-class of cryogenic labels. The DYMO LabelWriter 550 has 300 DPI Resolution, Includes Power Adapter, EMEA Power Cable, USB Cable, Quick Start Guide, Sample Roll and DYMO Connect for Desktop Software for Windows and Mac</p>	N/A	N/A	1 inch	1 printer	4 179 Kč
PKT-LWCS-505-EMEA	<p>Our DYMO® LabelWriter™ 550 printing kit provides an affordable, complete printing solution for labs on a budget. Supplied with the DYMO LabelWriter™ 550 direct thermal printer, free DYMO Connect label design/barcoding software, and two rolls of LWCS-class labels, this kit provides you with all you need to start printing labels for your lab.</p> <p>What's included:</p> <ul style="list-style-type: none"> - LabelWriter™ 550 Label Printer - Starter roll of address labels - Two rolls of LWCS-505-WH - Power Adapter - Power Cable - USB Cable - Quick Start Guide 	1" x 0.5" / 25.4mm x 12.7mm	N/A	1 inch	1 kit	6 423 Kč
PKT-LWCS-506-EMEA	<p>Our DYMO® LabelWriter™ 550 printing kit provides an affordable, complete printing solution for labs on a budget. Supplied with the DYMO LabelWriter™ 550 direct thermal printer, free DYMO Connect label design/barcoding software, and two rolls of LWCS-class labels, this kit provides you with all you need to start printing labels for your lab.</p> <p>What's included:</p> <ul style="list-style-type: none"> - LabelWriter™ 550 Label Printer - Starter roll of address labels - Two rolls of LWCS-506-WH - Power Adapter - Power Cable - USB Cable - Quick Start Guide 	0.94" x 0.5" + 0.375" / 23.9mm x 12.7mm + 9.5mm	N/A	1 inch	1 kit	6 423 Kč